

*Diagnostic • Surgical • Endoscopic Lamps
Lamps pour endoscopies • diagnostic et chirurgie*

MediLUME

For Zeiss model
304949-9001-000

Lamp
Order
Code:

ML-VAR700R

A Solution for the Zeiss OPMI VARIO 700 Surgical Microscopes

This MediLUME kit allows you to quickly and easily refurbish your 300 watt xenon lamp models to full OEM performance.

Each kit contains a pre-assembled heat-sink assembly with a new xenon lamp, a printed circuit board and step-by-step instructions.

To place an order or for more information please visit our website at:

www.medilume.com

MediLUME

Disclaimer: MediLume is not affiliated with Zeiss and use of any product names or codes are used solely to illustrate physical and technical compatibility of a given lamp or module with the respective equipment.